PHILIPS

Connected Blood Glucose Monitoring System

User manual

All images for reference only.

ART-303 Rev. 4

Icon definitions

Symbol	Description
Carrier Name	The name of the cellular carrier to which the meter is currently connected.
	Indicates that the current battery charge status. The battery energy remaining is depicted in a standard manner, familiar to most users, in which red indicates near-total discharge, which additional levels of white indicate progressively more complete charge states.
	Indicates that the battery is actively charging. These items are displayed in an animated fashion from top to bottom in order to depict active charging of the meter battery.
	Indicates the battery is getting low and needs to be charged.
	Indicates the strength of the wireless signal between base station and the meter in a manner familiar to most users. Additional bars of signal depict increased signal strength. Users who have two bars or less may experience slow and unreliable transmission and reception of data.
	The triangle icon indicates that the meter is roaming on the network of a cellular carrier different from the home carrier.
	The letter "X" is used to indicate no wireless signal connectivity.

This product contains a rechargeable internal battery. To conserve battery power this product will turn off the screen after one minute of non-use. To turn on the product either insert a test strip or press the Home or Power button.

Not intended for sale within the EU.

Contents

A. Intended use	4
B. Important safety instructions	4
C. Getting started	6
D. Your Connected Blood Glucose Meter (BGM)	8
D.1 Testing your blood glucose	12
D.2 Testing with control solution	18
D.3 Functions and settings	22
E. Customizing your meter	33
F. Classic mode	39
F.1 Testing your blood glucose	40
F.2 View your history on the meter	43
F.3 Testing with control solution	43
F.4 Setting up your BGM	47
G. Maintenance of your BGM and Lancing Device	48
H. Storing Your BGM and Test Strips	53
I. Troubleshooting	54
J. About your system	60
K. Terms of Use	65

A. Intended use

The Philips Connected Blood Glucose Monitoring System is intended for the quantitative measurement of glucose in fresh capillary whole blood samples drawn from the fingertips, forearm, or palm. It is intended for lay use by persons with diabetes to aid in diabetes management. It is indicated for use at home (over the counter [OTC]), and should be used only by a single patient, and should not be shared. Test results are plasma-calibrated for easy comparison to clinical laboratory methods. Testing is done outside the body (in vitro diagnostic use). The Connected Blood Glucose Monitoring System consists of the Philips Connected Blood Glucose Meter (BGM), Philips Test Strips, and Philips Control Solutions. The Connected Blood Glucose Monitoring System is not indicated for the diagnosis or screening of diabetes or for neonatal use. Palm and forearm testing should be done only during steady-state times when glucose is not changing rapidly. The BGM uses cellular data transmission to send test results to the Philips Virtual Care Portal, and to receive messages from the Virtual Care Portal. The Connected Blood Glucose Monitoring System is not intended to provide automated treatment guidance or decisions, nor is it to be used as a substitute for professional healthcare judgment.

B. Important safety instructions

The Connected Blood Glucose Monitoring System and lancing device are for single patient use.

- Do not share them with anyone, including family members!
- Do not use on multiple patients!
- Do not transfer ownership.

Because all parts of the kit may come in contact with your blood, all parts are considered biohazardous. Kit parts can potentially transmit infectious diseases, even after cleaning and disinfection. Dispose of in accordance with local regulations.

Users should wash hands thoroughly with soap and water, and dry thoroughly after handling the meter, test strips, and any lancing device.

For cleaning and disinfecting procedures to prevent the spread of diseases such as Hepatitis B, see "Cleaning and disinfecting your BGM" (pg. 50 ).

Use of this device on multiple patients may lead to transmission of Human Immunodeficiency Virus (HIV), Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), or other bloodborne pathogens.

**For technical support, contact Customer Care at
(888) 666-4530 (US) or (1) (978) 610-4688 (International)
Hours of operation: Monday - Friday, 9AM - 9PM ET.**

References

- "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)
<http://wayback.archive-it.org/7993/20170111013014/http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
- "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010)
<https://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>
- "Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition. Clinical and Laboratory Standards Institute (CLSI), document M29-A3, (ISBN 1-56238-567-4).

C. Getting started

Thank you for choosing the Connected Blood Glucose Monitoring System. The system is comprised of the Connected Blood Glucose Meter (BGM) and the My Virtual Care portal (when applicable). The meter provides accurate readings sent immediately to your online account to store, track, and chart. Personalized messages may also be sent to the meter. This manual will provide you with all you need to know to use your new meter. So let's get started!

Unboxing the BGM for the first time

Before we begin, take a moment to go over the contents of your BGM kit. Your kit should include:

Philips Connected
Blood Glucose Meter
(battery included)

AC Charger


Nylon Case


Lancing Device

(instructions for use included)

All images are for
reference only

User Manual

Quick Start Guide

The BGM kit does not contain testing supplies.

The BGM is compatible with Philips or BioTel Care® Test Strips and Philips Control Solution or BioTel Care® Glucose Control Solution only.

The following supplies are provided separately:

- Test Strips
- Lancets
- Control Solution (L1 and L2)

Using your new BGM:

- Make sure you charge your meter for at least 6 hours before first use.
- To charge your meter, use the AC charger included in your kit.

Important:

- Please remove the protective film over the screen.
- Wait for the meter to finish the boot process to start a test when it powers on.

You should use the BGM for:

- Testing blood glucose
- Testing only with fresh capillary whole blood samples

You should NOT use the BGM for:

- The diagnosis of Diabetes
- Testing blood glucose on newborn infants
- Testing of arterial blood
- Testing of venous blood

D. Your Connected Blood Glucose Meter (BGM)

1 POWER button ⏻

- Short-press (press and release within 2 seconds) of power button turns on/off the screen
- Long-press (press and hold for more than 2 seconds) of power button when screen is on pops up a window with the following options: Power Off, Restart, and Airplane Mode / Silent, Vibrate, and Sound
- Long-press of power button will turn on the meter when powered off

2 Up/down

- Buttons to scroll through menus

3 Enter button

4 Display touch screen

5 Test strip port (white side up)

6 Battery access notch


- Use the notch to remove the battery cover and replace the battery

7 Home button

- Press of home button returns to the home screen except in Test Mode

8 Charging port

Philips Test Strips

Important notes on test strips:

- Because of the auto coding feature of this meter, the strip insertion will be snug. This is to be expected and users find it easy after a few attempts.
- Check the expiration date printed on the test strip bottle and its package. DO NOT use expired test strips.
- Close the bottle immediately after taking out a test strip.
- Use each test strip immediately after removing it from the bottle.
- Do not use wet or damaged test strips.

The Home screen

On the Home screen, you can see the following:

- ① Wireless signal
- ② Battery level
- ③ Current time
- ④ Last test reading value
- ⑤ Date and time of last reading
- ⑥ Available widgets indicator
- ⑦ Applications Bar
- ⑧ Application Drawer

On the Home screen, swipe from the left edge of the screen to right or press the menu button in the top left, to open the Application Drawer to access most of the meter's options:

- Home
- View/Add Logs
- Messages
- App Store
- Settings
- Control
- Solution Test
- Test Mode
- Summary
- Applications
- Reminders
- Support

The Application Drawer has access to many of the apps available on your BGM. To access all available applications, press the Apps icon, or select using the UP/DOWN/ ENTER buttons, on the bottom Apps Bar of the Home screen.

Preparation and safety

Before you begin, make sure to:

- Set up your BGM properly.
- Review the Quick Start Tutorial that appears on screen upon the first time powering on the device. The tutorial shows how to use the BGM and several of its main features.

Note: The Quick Start Tutorial may be reopened again through the Support App.

- Run a control test before first use. (See "Step-by-step guide to control solution testing" on pg. 19 .)
- Wash your hands and testing site thoroughly with warm water and soap. Dry your hands well.
- Test in a place where the temperature range is between 50°F (10°C) and 104°F (40°C). (See "Troubleshooting" on pg. 54 .)
- Keep away from direct sunlight and heat. Store the test strip bottle in a dry cool place.

- Record the “date opened” on the bottle label. Discard the bottle and any remaining test strips 6 months after the date of opening, or if expired, which ever comes first.
- Do not operate the meter outside the temperature range of 50°F (10°C) to 104°F (40°C).

You will need:

- Your BGM
- A new Philips or BioTel Care® Test Strip
- A lancing device with a clean, unused lancet (or safety lancets)

IMPORTANT WARNINGS ABOUT INACCURATE RESULTS

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.

FOLLOW THESE IMPORTANT SAFETY GUIDELINES

- Please follow cleaning and disinfecting instructions for the meter and lancing device to prevent the risk of bloodborne pathogen transmission! (pg. 50 
- Follow any cleaning/disinfection instructions that may come with the lancing device you use.
- The meter and lancing device are intended for a single user and should not be shared.
- Make sure you wash with soap and water and dry hands thoroughly before and after handling the meter, lancing device, or test strips.

D.1 Testing your blood glucose

1. Power on the BGM with the **POWER** button on the right side of the BGM. Long press if it is turned off, short press if it is in SLEEP mode. The device will automatically display the Home screen.
2. Press on the TEST MODE icon or use the three buttons on the left side of the BGM to enter BG Test Mode.
When ready to check blood, "Insert a new test strip" will be displayed.

3. Take out a test strip from the bottle, and close the bottle immediately. (Some meters require the user to enter a code from the vial. Philips and BioTel Care® Test Strips do not require a code.)
4. Firmly insert the test strip into your BGM until it is fully seated. The Test Strip Port light will turn on.

5. Prepare your lancing device (refer to the lancing device's Instructions for Use).
6. Lance your finger (refer to the lancing device's Instructions for Use). Apply gentle pressure near the puncture site if needed, until a blood drop appears.
7. When finished, the lancet should be removed and disposed in accordance with local waste regulations.

Important information for the lancing device

- Use the lancing device included in your Connected Blood Glucose Monitoring System to obtain a blood sample. The minimum blood sample volume is 0.8 microliters.
- Purchase an additional lancing device by calling Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International), (9AM - 9PM ET, M-F).
- If using another type of lancing device, please refer to the instructions provided with the lancing device.

Important safety procedures regarding lancets

- To avoid contamination, use a new lancet every time you test. It is important to discard the used lancet carefully after each use to avoid unintended lancet stick injuries. Used lancets may be considered biohazardous waste in your area. Be sure to follow your local regulations for proper disposal.
- If alcohol wipes are used to clean your finger or palm, make sure your skin is dry before the blood sample is obtained.
- Lancing devices should NEVER be used by more than one person.
- If the seal is broken, discard the lancet.


How to choose the cap

- For fingertip testing, use the colored cap included on the lancing device.

- For alternate site testing on your palm, use the clear cap. This package does NOT include an alternate site cap. If you require one, please contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM - 9PM ET, M-F).

Tips for fingertip testing

- Gently massage your hand and finger toward the puncture site to form a drop of blood. Do not "milk" or squeeze around the puncture site.
- Lance the side of your fingertip to avoid soreness. To avoid calluses, choose a different lancing site each time. Position the end of the adjustable comfort tip against the side of the fingertip. Press the trigger button, and then lift the lancing device away from the finger after the puncture is complete.

Testing using an alternate test site

Important information on alternate site testing

Your BGM is capable of testing your blood glucose from your palm, usually referred to as an alternate site test, or AST.

Alternate site testing can be less painful than fingertip testing. However, because of the physiological differences between your fingertip and palm, alternate site testing results may be significantly different.

Alternate site testing should not be used to calibrate a continuous glucose monitor nor should it be used in insulin dose calculations.

Alternate site testing should only be used during steady state blood glucose conditions when glucose is not changing rapidly.

You should consult with your health care provider before using alternate site testing.

Testing on your palm

Prepare the lancing device with a new lancet and clear cap. Cock the lancing device so that it is ready to use. Select a site on your palm that does not have veins, moles, excessive hair, or bone. Vigorously rub the side of your palm with your other hand for at least 10 seconds. You should feel it get warm.

Press the opening of the clear cap against a fleshy area of the palm. Press the lancing device up and down in a slow pumping motion for several seconds to help the flow of blood.

Keep steady pressure on the area and press the release button. Do not be afraid to apply pressure to the area with the lancing device to assist the flow of blood. Continue to press the cap against the area for a few seconds to allow blood to come to the surface. You can see the blood drop through the clear cap. Lift the lancing device straight up. Do not smear the blood when removing the lancing device. If the blood drop is too small, the meter will show **Error: Not Enough Blood**. Insert a new strip and obtain a new blood drop.

If the results from alternate site testing differ from how you feel, use fingertip testing instead.

8. Applying blood to the test strip:

Gently touch the test strip to the drop of blood. The test strip will absorb the blood. Keep the test strip in the blood drop until you see the blood sample end of the test strip fill.

9. Tagging your test result:

The screen will start to count down. After 6 seconds, your glucose level will appear on the display screen. You can select a Reading Tag (e.g., "Before Breakfast") from the list displayed during or after countdown. Swipe down or use the **DOWN** key to select the correct tag. Press on the **CONTINUE** button or use the **ENTER** key to confirm. Your test result will be automatically stored into the meter history. If a Reading Tag is not selected in the allotted time, it will automatically be tagged as "No Selection".

10. Sending your results:

After pressing the **CONTINUE** button, your BGM will send the results to the portal.

11. Remove the used test strip and discard it properly.

Clinical response screen:

During data transmission, you will see a display that compares your reading to your set target. To change your set target, call your health care provider (if already connected) or Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F). You may also see helpful guidance from your health care provider regarding this reading to manage your diabetes.

Important steps to follow if you see High or Low

If you see High or Low displayed, your blood glucose level may be above 600 mg/dL or below 20 mg/dL. Wash your hands. Dry your hands completely. Test again using fingertip testing. DO NOT test on your palm or forearm. If you still receive the same message, test your BGM with a control solution. Refer to "Step-by-step guide to control testing." (pg. 19 ).

Optional: After successful data transmissions to the portal, you may receive a multiple choice question to answer. This may be a question from your health care provider or relate to reordering supplies. Swipe up/down or use the **UP/DOWN** buttons to select the appropriate answer. Press on the option or use the **ENTER** button to transmit.

If your data transmission is not successful

If the device cannot communicate with the portal, the following message will appear:

No network connection. Reading will be sent when a connection returns.

Your test results will be stored on your device. You can send the reading on your own by going into Logs from the main menu. Select the unsent reading displayed in Logs and press the Sync button to transmit. For the transmission to be successful, network connection is required.

The meter will attempt to forward unsent readings automatically every 15 minutes, or following the next sent reading. Readings that have not been transmitted to the portal will be marked with a yellow circular hazard icon in the Logs list.

IF DATA TRANSFER FAILS, YOU MAY USE THIS PRODUCT LIKE ANY OTHER BLOOD GLUCOSE METER. FOLLOW YOUR HEALTH CARE PROFESSIONAL'S INSTRUCTIONS ON TREATMENT.

Understanding your blood glucose test results

Blood glucose value will vary depending on food intake, medication, health, stress, and exercise. The ideal ranges for people without diabetes are as follows¹:

- 70 to 100 mg/dL before meals
- 70 to 140 mg/dL after meals

What to do if you get high or low readings AND you feel ill.

If the BGM displays results that are "High" or "Low," or you get a result that is more than 250 mg/dL or below 50 mg/dL AND you feel ill:

- Treat your diabetes according to the instructions from your health care provider.

If the BGM displays results that are "High" or "Low," or you get a result that is more than 250 mg/dL or below 50 mg/dL AND you DO NOT feel ill:

- Test your BGM with a control solution: refer to "Step-by-step guide to control testing." (pg. 19 )
- Wash your hands, dry thoroughly and test again using your fingertip.

If you still get a high or low reading, contact your health care provider immediately.

References

1. American Diabetes Association. Standards of Medical Care in Diabetes—2023. *Diabetes Care*. 2023; 46:S1-S285.

D.2 Testing with control solution

Control solution testing will confirm that your BGM and test strips are working properly. You should perform control solution testing when:

- Using the meter for the first time
- Using a new bottle of Philips or BioTel Care® Test Strips
- You left the test strip bottle cap open
- You dropped the meter
- You suspect the meter and test strip are not working properly
- The test results do not reflect how you feel
- You want to practice taking a blood glucose test
- The test strips are exposed to environmental conditions outside the temperature requirements printed on the label

Purchase additional control solution by calling Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

Step-by-step guide to control solution testing

1. Power on the BGM with the **POWER** button on the right side of the BGM.
The device will automatically display the Home screen with the last reading as shown.

2. Press on the Control Solution Test icon or use the three buttons on the left side of the BGM to enter Control Solution Test Mode.

Switch mode to L1 Control Solution Mode.

To test level 1 (L1) control solution, press on the **L1 CONTROL** button or use the **UP/DOWN** buttons to highlight **L1 CONTROL** button.

Press on the Control Solution Test icon or use the three buttons on the left side of the BGM to enter Control Solution Test Mode.

Switch mode to L2 Control Solution Mode.

To test level 2 (L2) control solution, press on the **L2 CONTROL** button or use the **UP/DOWN** buttons to highlight **L2 CONTROL** button.

- 
- Take out a Philips or BioTel Care® Test Strip from the test strip bottle. Close the bottle immediately.

- 4. IMPORTANT:**

- Firmly insert the test strip into your BGM until it is fully seated. The Test Strip Port light will turn on. Continue to push the strip in until **"Apply Level 1 (2) Control Solution"** appears.
- The term "Philips" or "BioTel" should be facing up.
- The sample tip should be facing up.
- Do not touch the sample tip until applying control solution.

- Use level 1 solution for L1 mode. Use level 2 solution for L2 mode. Squeeze a drop of control solution onto a clean, dry, non-absorbent surface. Immediately replace the bottle cap on the control solution bottle. Do not apply control solution to the test strip directly from the bottle.

* Control Solution Range
(Images are for reference only. Actual numbers may vary.)

- Gently touch the tip of the test strip to the drop of control solution. You should see the sample end of the test strip fill.
- The screen will start to count down. After 6 seconds, the control test result will appear on the meter's display screen.
- Compare the reading on the screen to the "L1 Control Solution" or "L2 Control Solution" range printed on the test strip bottle or package. If the reading on the meter falls within the appropriate range, select the YES tag. If not, select the NO tag. If the reading falls outside the level printed on the bottle or package that you used, please see "Control solution troubleshooting" (see section below).

Control solution readings are stored on the meter but are not included in the Blood Glucose testing logs.

Important tips about control solution

1. Only use Philips Control Solution or BioTel Care® Glucose Control Solution with the Philips or BioTel Care® Test Strips. Other brands of control solution or test strips may produce inaccurate results.
2. Always check the expiration date. DO NOT use expired control solution.
3. DO NOT FREEZE. Store the control solution at room temperature, between 39°F (4°C) and 86°F (30°C).
4. Never touch the tip of the control solution bottle to the test strip. It will damage all remaining solution in the bottle.
5. Record the “date opened” on the bottle label.
6. Do not use the control solution if more than 3 months since first opened.

Control solution troubleshooting

If your control solution testing is out of range (too high or too low), it may be caused by the following:

- Wrong brand of control solution
- Control solution not at room temperature
- Expired or contaminated control solution
- Damaged test strip
- Touched the tip of solution bottle to test strip
- Meter malfunction
- Testing outside of meter operating temperatures (50°F – 104°F) (10°C – 40°C)

Check the expiration date and open bottle date on the control solution and on the test strips. Repeat testing using a new test strip. If the result is still out of range, use a new bottle of control solution and retest.

If the problem persists, call Customer Care at
(888) 666-4530 (US) or (1) (978) 610-4688 (International)
(9AM-9PM ET, M-F).

Do not perform any blood glucose tests until you resolve this issue.

D.3 Functions and settings

View your logs on the meter

View your logs on your meter

Your BGM can store at least 10,000 logs. Each log is stored with the date and time.

1. Power on the BGM by using the **POWER** button.
2. Press on the **LOGS** icon or use the **UP/DOWN/ENTER** buttons to access logs.
3. Swipe up/down or use the **UP/DOWN** buttons to scroll through stored readings and other entries. (Note: select any log to view more details)
4. Press the **Filter** button to select a subset of log types for display.
5. To exit the Logs menu, press the **HOME** button.

Blood glucose readings

Blood glucose readings display in a colored circle. Each color indicates in which target range the reading falls:

Blue: below target

Green: on target

Orange: above target

Red: high above target

Control solution and other Logs

- Control solution logs are displayed in an uncolored circle.
- Non-Glucose logs are displayed in a uncolored rectangular box.
- Blood Pressure logs from a connected device are displayed in a colored (green, yellow, orange, or red) circle.
- Manually entered Blood Pressure logs are displayed in a colored (green, yellow, orange, or red) square.

Control
Solution Log

Non-
Glucose Log

Blood
Pressure Log

Manually
Entered Blood
Pressure Log

List of Log icons

The following 14 different categories can be displayed in the Logs Screen: Blood Glucose, Carbs, Medicine, Activity, Weight, Blood Pressure, A1C, SPO₂, Temperature, Sleep, Body Mass Index, Water Intake, Heart Rate, and Meal. Each category has its unique icon as shown in the picture.

● Blood Glucose	● SPO ₂
† Carbs	● Temperature
● Medicine	● Sleep
● Activity	● Body Mass Index
● Weight	● Water Intake
● Blood Pressure	● Heart Rate
● A1C	● Meal

Symbol	Descriptions
	The keyboard icon is used to indicate that a log has been manually entered.
	The blood drop icon is used to indicate that a Blood Glucose log has been automatically added through BG Test Mode.
	The blood pressure cuff icon is used to indicate that a Blood Pressure log has been automatically added from a connected device.
	The weight scale icon is used to indicate that a Weight log has been automatically added from a connected device.
● ● ●	This icon is used to display the availability of a drop-down menu
	The thermometer icon is used to indicate a Body Temperature log.
	The pulse oximeter icon is used to indicate an SPO ₂ log.
	This icon indicates that the associated log has not been uploaded to the portal. You can select the log and retry to send it.

Add Logs on the meter

Add your logs on your meter

Your BGM allows you to manually add Blood Glucose, Carbs, Medicine, Activity, Weight, Blood Pressure, A1C, SPO₂, Temperature, Sleep, Body Mass Index, Water Intake, Heart Rate, and Meal.

1. You can enter the **Add Log** mode from the Logs Screen or Applications Menu by pressing on the **Add Log** button.
2. Select the log type by pressing on the option or use the **UP/DOWN** button.
3. Press the **CREATE A LOG** button or use the **ENTER** button to enter the Add Log Screen where you can enter log details.
Press the **EDIT** button to edit or create a new log type.
(See step 5)
4. You can enter the log details such as amount, tag and timestamp and press **SUBMIT** to save and upload this log to the portal.

Create new Log Type

1. You can enter the Edit Logs screen from the Add Log screen by pressing the **EDIT** button. Select a Log Type (Carbs, Medicine, or Activity) and press the **NEW LOG TYPE** button, where you can enter the new type name and unit. Upon entering the new type name and unit, press **SAVE**.

Examples of new Log Types: Carbs, Breakfast, Lunch, Dinner, Medicine, Insulin, Vitamins, Activity (walking, jogging, biking, etc.)

2. Press **SAVE** on the Add Log screen to save the log type.
3. Delete new Log Type:
Once a new Log Type is created, it may be deleted from the List of Logs. From the Add Log screen, select the new Log Type to be deleted. Once highlighted, press the **EDIT** button. When the new Log Type screen appears for editing, press the Delete button. When the Add Log screen appears again, press the **SAVE** button.

View messages

Display messages on your meter

Your BGM allows you to check historical messages that have been received after each blood glucose test. The messages will be displayed in chronological order with most recent at the top.

1. You can enter the Messages Screen from the Applications Menu and pressing the Messages icon or use the buttons to select and enter. You may also use the Application Drawer to select Messages.
2. You can swipe up/down or use the **UP/DOWN** buttons to navigate through each message.
3. To return to the Home screen, press the **HOME** button.

Summary of past readings

Display statistics of your past readings on your meter

Your BGM allows you to check the statistic result of your past blood glucose readings. The summary screen has 5 graphs: Personal Testing Goal, Readings by Range, Hourly Averages, Readings out of Range, and Average Result by Tag.

1. You can enter the Summary screen from the Applications Menu by pressing the Summary icon or use the buttons to select and enter. You may also use the Application Drawer to select Summary.
2. You can swipe left/right, press the "<" and ">" icons or use the **UP/DOWN** buttons to navigate to each graph.
3. You can also select one of the following timeframes to generate the graph: Today, 7, 14, 30, 90 days, All.
4. You can view the graph legend by selecting the **LEGEND** button.
5. To return to the Home screen, press the **HOME** button.

Reminders

Your BGM allows you to create and manage reminders for blood glucose readings, medications or activities, etc.

1. You can enter the Reminders screen from the Applications Menu by pressing the **REMINDER** icon or use the buttons to select and enter. You may also use the Application Drawer to select Reminders.

2. You can create a new reminder by pressing the **CREATE REMINDER** button.

You can edit an existing reminder by pressing on it; this will bring up the reminder detailed information.

In the detailed information, you can change the time, title, and days which the reminder are active by pressing on each area. Pressing the **SAVE** button saves the reminder setting. Pressing the **DELETE REMINDER** button deletes the reminder.

3. To return to the Home screen, press the **HOME** button.

Support

Get on-screen support from tutorials on your meter.

Your BGM allows you to view video tutorials, the Quick Start Tutorial, and the user manual, all on your meter.

1. You can enter the Support screen from the Applications Menu by pressing the Support icon or use the buttons to select and enter. You may also use the Application Drawer to select Support.

2. You can open a page-by-page slideshow of the user manual by pressing the **VIEW** button or use the **UP/DOWN** buttons to highlight user manual and then press the **ENTER** button. You can swipe up/down and press on the tutorial you want to view or use the **UP/DOWN** buttons to highlight the tutorial and then press the **ENTER** button.
3. To return to the Home screen, press the **HOME** button.

Setting up your BGM

Your BGM time is automatically updated with the local date and time when a cellular connection is established.

If you would like to view your meter's preset settings, navigate to the Settings screen.

Navigating to the Settings screen

1. Power your BGM on.
2. You can enter the Settings screen from the Applications Menu by pressing the Settings icon or use the buttons to select and enter. You may also use the Application Drawer to select Settings.
3. Press on the **SETTINGS** icon or use the **ENTER** button to enter the Settings Menu.

Setting options

Volume

Change the volume of your meter by moving the slider thumb.

Brightness

Change the brightness of your meter's screen.

Airplane Mode

Turning on Airplane Mode disables cellular network access. You will not be able to transmit your readings while Airplane Mode is turned on.

Your meter's clock will continue to work properly while in this mode.

If a time zone change occurs, your meter will not automatically update to the new time until Airplane Mode has been disabled. The meter will reset the time from the new network.

Languages

Change the language displayed by your meter.

Time format

Change the time display as 24-hour or 12-hour.

Date format

Change the date display format.

Tag timeout

BG tag timeout determines the duration of the tagging time before a new blood glucose reading is automatically sent to the portal if no selection is made.

Reminders

Press here to manage test reminders.

Personal Testing Goal

Display and set your goal for number of tests per day.

Mode

This changes the user interface styles among "Normal" and "Classic" views.

Theme

This changes the theme (color) of the user interface. The three options are blue, white, and gradient.

Blood Glucose display

This changes the “Blood Glucose Last Result” widget display on the Home screen. The three options are horizontal, vertical, and graphical.

Meter Name

You may enter a name on this screen that can help identify your meter and will be displayed on the top of your Application Drawer and Home screen.

Check connectivity

You may check whether it connects to the portal.

About

This screen provides a listing of device software and related information about your meter. This is helpful if you need technical support.

Version information

Press here to check the current application versions.

Apps

Download software updates for your BGM from the App Store.

Your BGM allows you to download software updates available for your BGM from the App Store. You can enter the App Store screen from the Applications Menu by pressing the App Store icon or use the buttons to select and enter. You may also use the Application Drawer to select App Store.

Critical updates

When a critical software update is available, an **INSTALL NOW** prompt will appear on your screen. This prompt will take you to the App Store to select and install the update. You will not be able to use other meter functions until a critical update has been installed.

Non-critical updates

When a non-critical software update is available, an **INSTALL NOW** or **TRY LATER** prompt will appear on your screen. The **INSTALL NOW** prompt will take you to the App Store to select and install the update. The **TRY LATER** prompt will display the software update prompt again after 24 hours, and continue every 24 hours until the update has been completed.

E. Customizing your meter

Home screen Application Bar customization

Applications may be chosen to be displayed on the Applications Bar of the Home screen. This option allows for improved ease-of-use and quick access to favorite applications.

1. Press the Apps Menu button on the bottom of the Home screen to navigate to the Applications page (refer to Applications Menu image).
2. Press the three dots in the top right corner of the applications screen for the application customization drop-down menu (refer to Customization Menu image).

Home screen

Applications Menu

3. To customize the Home screen Applications Bar, select the Home screen Apps button from the drop-down menu (refer to Home screen Apps image).

4. On the Select Home screen Apps page, change up to 3 additional available applications to display on the Home screen App Bar. No more than 5 total applications may be selected at any time - Test Mode and Applications Menu icons are always selected by default. Select the Done button to Set Applications (refer to Home screen Apps image).

Widget customization


For quick access to application functions through mini-applications, widgets are made available to be used on the Home screen by swiping Left/Right. Select applications support widgets and may be added to the Home screen.

1. To customize the Home screen Widgets, select the **Home Screen Widgets** button from the drop-down menu (refer to Customization Menu image).
2. On the Create Home screen Widget page, select available applications to display on the Home screen. Select **Done**

button to choose Widgets. Drag Apps into desired order for display and select **Done** (refer to Home Screen Widgets and Widget Order images).

Home Screen Widgets

Widget Order

Delete applications

With the addition of new applications, user-selected applications may be allowed to be removed based on use or personal preference.

1. To remove applications, select the **Delete Apps** button from the drop-down menu (refer to Customization Menu image).
2. On the Delete Selected Apps page, if available for removal, press the "X" in the top right corner of the app icon to delete from the Applications screen. Select **Done**.

Note: Deleted applications may reappear following system updates (refer to Delete Applications image).

Delete Applications

Blood Glucose Last Reading display customization

To better understand the Last Reading taken on the meter, a visual representation of the ranges is displayed. The Blood Glucose Display is able to be customized for improved understanding and clarity.

1. Select the Settings application either via the Application Drawer or Applications Menu (refer to Settings Icon image).
2. Scroll down and select the Blood Glucose Display button.
3. From the drop-down menu, select one of the following Last Reading Gauges: Graphical (default), Vertical, or Horizontal.

Settings Icon

Graphical

Vertical

Horizontal

Theme customization

For better viewing and personal preference, the theme may be customized. The varying themes change the background color and style of the meter.

1. Select the Settings application either via the Application Drawer or Applications Menu (refer to Settings Icon image).
2. Scroll down and select the Theme button.
3. From the drop-down menu, select one of the following Themes: Blue, White, or Gradient (default).

Other apps

Along with the BGM default apps, there are other apps that may appear on your device. These apps include:

Connected Devices

You can connect select devices (ex: pulse oximeter, blood pressure monitor or weight scale) to your BGM via Bluetooth and send the results right to your meter. After choosing the Connected Devices app icon on your meter, select the Scan button in the top left of the screen to pair your device (if this is your first time in the app, you will be asked to turn on Bluetooth via a pop-up, click Yes). Follow instructions sent with your connected device to pair the device with your meter. Once paired, you may take a test with your connected device. Connected device results are synced with the app and displayed on your BGM screen. The result is sent to the portal and saved as a log. If needed, user manuals may be accessed by pressing the Info button in the top right of the main Connected Devices screen.

Survey App

To better engage with you, surveys may be sent to your meter to gain your feedback. Upon receipt of a new survey, a pop-up will appear on the screen. You can choose to either open the survey or dismiss the pop-up. If you dismiss the survey, it will still appear in the Survey App. To access the new survey, you may either select Open on the pop-up or go to the Survey App by pressing the Survey icon on your meter.

In the app, you can select any new or listed surveys. Once in a survey, you may answer the on-screen questions. After you have answered all survey questions, your meter will submit the results to the portal. Once a finished survey is sent to the portal, the results will not be stored on your meter.

Note: If a survey is closed before it is finished, your answers will not be saved.

Alert App

The Alert App is used to send a direct message to your device from your program or health care provider, in order to provide you with feedback or questions. Alerts are sent to your meter and are stored within the Alert App. Once an alert is received, a pop-up will appear on your device that may be read and dismissed. If a response is needed for an alert, it will be noted in the Alert App and you can select an answer to send. Prior Alerts are stored in the app after being read and dismissed.

Note: Not all Alerts require a response.

Refill App

When you are running low on testing supplies for your BGM, you may be able to request refills right from your meter. Once opened, the Refill App will display on-screen prompts and questions to complete a Refill request. After the questions are answered, your request will be sent to our system for review and order placement. Once your request is received, you may receive a follow up call if there are any additional questions.

Note: If you exit the Refill request before you answer all questions, the answers will not be saved or sent to Philips.

Note: If an app appears on your BGM that is not listed here, please contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F) or VCM-customercare@philips.com for more details. Not all apps listed in this section are offered to all users.

F. Classic mode

You can select the Classic mode under the Mode option on the Settings screen, and the Classic style Main Menu screen will be displayed. It mimics the User Interface in previous models.

Main Menu

In this Home screen, swipe up/down or use the **UP/DOWN** buttons (located on the left of your meter) to access the meter's options:

- Test Mode
- L1 Control Solution Mode
- L2 Control Solution Mode
- History
- Settings Menu


F.1 Testing your blood glucose

1. Power on the BGM with the **POWER** button on right side of meter. The device will automatically start in the Main Menu. Select Test Mode to test your blood glucose. When ready to check blood, "Insert test strip" will be displayed.

On the Test Mode screen, you can see the following:

- ① Wireless signal strength
- ② Battery level
- ③ Current time
- ④ Date and time of last reading
- ⑤ Last test reading

2. Take out a Philips or BioTel Care® Test Strip from the bottle, and close the bottle immediately. (Some meters require the user to enter a code from the vial. Philips and BioTel Care® Test Strips do not require a code.)
3. Firmly insert the test strip into your BGM until it is fully seated.
4. Prepare your lancing device (pg. 13 ).

5. Fingertip testing*:

Position the end of the adjustable comfort tip against the side of the fingertip. Press the trigger button, and then lift the lancing device away from the finger after the puncture is complete.

6. Testing using an alternate test site*

7. Applying blood to the test strip:

Gently touch the test strip to the drop of blood. The test strip will absorb the blood. Keep the test strip in the blood drop until you see the blood sample end of the test strip fill.

8. Tagging your data:

The screen will start to count down. After 6 seconds, your glucose level will appear on the display screen.

You can select a Reading Tag (e.g., "Before Breakfast") from the list displayed during or after countdown. Use the **DOWN** button to scroll to the correct tag. Press the **ENTER** button to select that tag.

Your test result will be automatically stored into the meter history. If a Reading Tag is not selected in the allotted time, it will automatically be tagged as "No tag".

9. Sending your results:

As soon as you select the tag, your BGM will send the results to the portal.

10. You can remove the test strip. Discard the used test strip properly.

* Please see details in Section D.1 - Testing your Blood Glucose (pg. 12)

Important steps to follow if you see High or Low

If you see High or Low displayed, your blood glucose level may be above 600 mg/dL or below 20 mg/dL. Wash your hands. Dry your hands completely. Test again using fingertip testing. DO NOT test on your palm or forearm. If you still receive the same message, test your BGM with a control solution: refer to "Step-by-step guide to control solution testing" (pg. 44 ).

OPTIONAL: After successful data transmissions to the portal, you may receive a multiple choice question to answer. This may be a question from your health care provider or relate to reordering supplies. Swipe up/down or use the **UP/DOWN** buttons to select the appropriate answer. Press on the option or use the **ENTER** button to transmit.

Clinical response screen

After data transmission, you will see a graph that compares your reading to your set target. To change your set target, call your health care provider or Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F). You may also seek helpful guidance from your health care provider regarding this reading to manage your diabetes.

F.2 View your history on the meter

View your logs on your meter

Your BGM can store at least 10,000 logs. Each result is stored with date and time.

1. Power on the BGM by using the **POWER** button.
2. Use the **UP/DOWN** keys to access History.
3. Use the **UP/DOWN** keys to scroll through the stored readings.

To exit the History menu, scroll to the bottom of your reading history and select “**BACK**”.

F.3 Testing with control solution

Control solution testing will confirm that your BGM and test strips are working properly. You should perform control solution testing when:

- Using the meter for the first time
- Using a new bottle of test strips
- You left the test strip bottle cap open
- You dropped the meter
- You suspect the meter and test strip are not working properly
- The test results do not reflect how you feel
- Practice taking a blood glucose test
- The test strips are exposed to environmental conditions outside the temperature requirements printed on the label.

Purchase additional Control Solution by calling Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM EST, M-F).

Step-by-step guide to control solution testing

1. Power on the BGM by using the **POWER** button on the right side of the BGM.

The device starts on the Main Menu. Once start up is completed, press Test Mode and the device will be ready to check blood glucose. You will see the "Insert Test Strip" screen as shown.

2. To access the main menu, use the **UP/DOWN** keys on the left side of the BGM. Switch modes from Normal Test Mode to L1 Control Solution Mode.

To test level 1 (L1) control solution, use the **UP/DOWN** buttons to highlight **L1 mode**. Press the middle (ENTER) button to select L1 mode.

Upon successful completion of level 1 (L1), test level 2 (L2) control solution:

Use the **UP/DOWN** keys on the left side of the BGM to access the main menu. Switch modes from Normal Test Mode to L2 Control Solution Mode.

To test Level 2 (L2) control solution, use the **UP/DOWN** buttons to highlight **L2 mode**. Press the middle (ENTER) button to select L2 mode.

3. You can select any mode you want by pushing the **ENTER** button.
4. Take out a test strip from the test strip bottle. Close the bottle immediately.
5. **IMPORTANT:**
 - Firmly insert the test strip into your BGM until it is fully seated.
 - The term “Philips” or “BioTel” should be facing up.
 - The sample tip should be facing up.
 - Continue to push the strip in until **“Apply Level 1 (2) Control Solution”** appears.
6. Do not touch the sample tip. Use Level 1 solution for L1 mode. Use Level 2 solution for L2 mode. Squeeze a drop of control solution onto a clean, dry, non-absorbant surface. Immediately replace the bottle cap on the control solution bottle. Do not apply control solution to the test strip directly from the bottle.
7. Gently touch the tip of the test strip to the drop of control solution. You should see the sample end of the test strip fill.
8. The screen will start to count down. After 6 seconds, the control test result will appear on the meter’s display screen.
9. Compare the reading on the screen to the “L1 Control Solution” or “L2 Control Solution” range printed on the test strip bottle or package. If the reading on the meter falls within the appropriate range, select the YES tag. If not, select the NO tag. If the reading

falls outside the level printed on the bottle or package that you used, please see "Control Solution Troubleshooting" (pg. 21 ).

Control Solution readings are stored on the meter but are not transmitted and displayed on the portal.

Important tips about control solution

1. Only use Philips Control Solution with Philips or BioTel Care® Test Strips. The BGM is compatible with the Philips Control Solution and BioTel Care® Glucose Control Solution, and Philips and BioTel Care® Test Strips. Other brands of control solution or test strips may produce inaccurate results.
2. Always check the expiration date. DO NOT use expired control solution.
3. DO NOT FREEZE. Store the control solution at room temperature, between 39°F (4°C) to 86°F (30°C).
4. Never touch the tip of the control solution bottle to the test strip. It will damage all remaining solution in the bottle.
5. Record the "date opened" on the bottle label.
6. Do not use the control solution if more than 3 months since first opened.

Control Solution Range (Images are for reference only. Actual numbers may vary.)

F.4 Setting up your BGM

Your BGM time is automatically updated with the local date and time when a cellular connection is established.

If you would like to view your meter's preset settings, navigate to the Settings screen.

Navigating to the Settings screen

1. **POWER** your BGM on.
2. Use the **UP/DOWN** buttons to reach the Settings Menu.

Setting Options

Alarms

Create up to 3 alarms to help remind yourself when to test next.

Date and Time

This page allows you to view the Date and Time settings of your device.

Volume

This page allows you to select the volume of your meter. Use the **UP/DOWN** keys to select the appropriate volume.

Brightness

Change the brightness of your meter's screen.

Languages

Change the language displayed by your meter.

BG Tag Timeout

BG Tag Timeout determines the duration of the tagging time before a new blood glucose reading is automatically sent to the portal if no selection is made (Note: Default BG Tag Timeout is 60 seconds).

Mode

This changes the user interface styles among "Normal" and "Classic" views.

Airplane Mode

Turning on Airplane Mode disables cellular network access. You will not be able to transmit your readings while Airplane Mode is turned on. Your meter's clock will continue to work properly while in this mode.

If a time zone change occurs, your meter will not automatically update to the new time until Airplane Mode has been disabled. The meter will update the time from the new network.

Personalize

You may enter a Name on this screen that can help identify your meter and will be displayed on your screen in Test, L1 and L2 mode. Use the **UP** and **DOWN** keys to choose a character and the **ENTER** button to enter it. When you are done entering your Name, press the **ENTER** button again after entering the last character.

About

This screen provides a listing of current software and related information about your meter. This is helpful if you need technical support.

G. Maintenance of your BGM and Lancing Device

Charging your BGM

Your BGM uses a rechargeable battery. It can be charged and discharged hundreds of times, but eventually the battery will lose its charge. If you notice significant reduction in the duration of your meter's operating time, please contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F), or e-mail VCM-customercare@philips.com.

If left unused, a fully charged battery will lose its charge over time. Always charge indoors. Under normal use conditions with good cellular signal strength (more than 2 bars), a fully charged battery should last up to 3 days.

To Charge Your Device

1. Plug the power cord into an electrical outlet.
2. Plug the other end of the cord into the BGM's charging port.

Charging Troubleshooting

If your device does not turn on when you press the power button on, your battery may not be fully charged. If after 6 hours of charging the device the device does not turn on, contact Customer Care at (888) 666-4530 or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

Warning:

- Please only use the included charger to charge your meter; other chargers may result in excessive heating of the device.
- Please call Customer Care and discontinue use if you feel the meter is heating up significantly when charging.

Cleaning and Disinfecting

What is the difference between cleaning and disinfecting?

Cleaning and disinfecting are not the same thing. Cleaning simply removes blood and dirt. Disinfecting kills germs and provides an extra level of safety. In the following procedure, the first disinfectant towel is used for cleaning, and the second towel provides for disinfection.

You should disinfect areas where there are both high concentrations of germs and a possibility that they will be spread to others. The Caltech DISPATCH Cleaner/Disinfectant with Bleach Premoistened Towels® have ingredients that clean and disinfect your device. While surfaces may look clean, many infectious germs may be lurking around. Some germs can live on surfaces for hours and even for days! The disinfection process shown in the pages that follow is effective in destroying dangerous Hepatitis B viruses.

If the device is being operated by a second person who is providing testing assistance to the user, the meter should be disinfected prior to assistance from the second person. See Cleaning your Lancing Device section (pg. 52 ) for instructions on cleaning the lancing device included with your kit. If using a different lancing device, see your device's user manual for decontamination/cleaning instructions.

To acquire canisters of Caltech DISPATCH towels, please contact:

Amazon - www.amazon.com

Sears - www.sears.com - (800) 697-3277

Cleaning and disinfecting your meter

Clean and disinfect your meter at least 1 time per week using Caltech DISPATCH Cleaner/Disinfectant with Bleach Premoistened Towels® using the following method:

1. Dispense the DISPATCH towel from its package or container.
2. Clean the outside of the meter to remove any blood or dirt.
3. Discard the used towel in a sealed container where it will not be touched by others.
4. Dispense a second DISPATCH towel from its package or container.
5. Disinfect the meter with the second towel by wiping all outside surfaces.
6. Discard the used towel in a sealed container where it will not be touched by others.
7. Wait at least one minute and make sure device is completely air dried before your next blood test.

The BGM has been validated for 1,825 cleaning and disinfection cycles. Please obtain a new meter after cleaning and disinfecting the meter 1,825 times or once the lifetime of the meter (5 years) has been reached, whichever comes first.

Meter cleaning do's and don'ts

- Do not allow cleaning solution to run in or around the test strip port. Doing so may cause a malfunction.
- Do not let the meter come in contact with water.
- Do not use any glass or household cleaners on the meter.
- **IMPORTANT:** Avoid getting cleaning solution, dust, dirt, blood, control solution, water, or any other substance in the meter's openings, including the test strip port, USB port and the test strip port LED light.
- **IMPORTANT:** If the meter is being operated by a second person who is providing testing assistance to the user, the meter should ALWAYS be disinfected prior to assistance from the second person.
- Stop using the meter and contact Customer Care immediately if you notice any of the following signs of deterioration:

- The meter does not turn on.
- The meter is damaged or has stopped working.
- The meter's screen has been cracked or has become clouded.

For technical support, contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

Cleaning your lancing device

Important: Do not immerse the body of the lancing device in water or any other liquid. Use isopropyl alcohol or soap and water to wipe the outside of the lancing device. If you wish, remove the cap, wash it in warm water, and rinse. Store lancing device in a carrying case or in a clean and dry environment.

IMPORTANT LANCING SAFETY MEASURES

If the lancing device is being operated by a second person who is providing testing assistance to the user, the lancing device and endcap should ALWAYS be disinfected prior to assistance from the second person.

Stop using the lancing device and contact Customer Care immediately for a replacement if you notice any of the following signs of deterioration:

- Any buttons have fallen off or stopped working.
- The device is cracked.
- The device is discolored.
- The device does not yield an adequate blood drop.

Warranty information

Please refer to the Terms of Warranty section (pg. 57  of this user manual. Please read this warranty in full. Your rights under the laws of your country are not affected.

H. Storing Your BGM and Test Strips

- Your meter is a delicate electronic device that should be handled with care. Severe shock, such as dropping the meter, could damage it.
- The meter and test strips only work in the temperature range of 50°F (10°C) to 104°F (40°C).
- Store the test strips between 39°F (4°C) and 86°F (30°C). Do not freeze. See the test strip instructions for more details.
- Do not leave the meter in a place that is extremely hot or extremely cold (e.g., anywhere near a heat source or in an extremely hot or cold car).
- Do not expose the meter or test strips to high levels of humidity, such as in the bathroom or kitchen.
- Immediately close the bottle cap tightly after removing a test strip.
- Do not attempt to open the meter. Opening the meter will void the warranty (unless replacing the battery or otherwise instructed by Customer Care) and may cause it to operate improperly. If you have technical problems or questions, turn your meter off and contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International). Support is available Monday through Friday, 9AM - 9PM ET.
- Dispose of in accordance with local waste regulations.

I. Troubleshooting

Error codes

In the event your device experiences an error, you will see the following screen along with an error message.

1. Error message: "Meter Failure"

What it means: As with any computer device, sometimes the meter may need to be restarted.

Action to take: Restart the meter and contact Customer Care if the problem continues. Call Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

2. Error message: "Used Strip Detected! Please insert a new strip."

What it means: Your test strip may be used or damaged.

Action to take: Repeat the test with a new test strip. Remember, add the sample only when you see the prompt on the display screen.

3. Error message: "Not Enough Blood" or "Not Enough Control Solution"

What it means: There is not enough blood or too much blood or control solution on the test strip to give a reliable reading.

Action to take: Remove the test strip and repeat the test with a new test strip.

4. Error message: "The temperature of the environment is too high!" or "The temperature of the environment is too low!"

What it means: The meter is too hot or too cold.

Action to take: Move the meter into an area that is between 50°F (10°C) and 104°F (40°C). After receiving this error, allow the meter to be at room temperature for at least 30 minutes before testing blood or control solution.

5. Error message: "No network connection. Reading will be sent when a connection returns."

What it means: The meter is unable to connect to the portal.

Action to take: See "Communication errors" below (pg. 56 ). If this is not resolved and you have more than 1 bar of cellular signal strength, call Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

6. Error message: "Reading is below the minimum 20 mg/dL. It may be caused by meter failure. Please insert another test strip and take a new reading."

What it means: The reading is below the minimum test range.

Action to take: Remove the test strip, insert another one and take a new reading. If the problem continues, please contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

Other messages

1. Icon message: Low Battery icon displayed

What it means: The meter is about to lose power.

Action to take: Charge the battery.

2. Icon note: Wireless signal indicator bars

Two bars or less of wireless signal strength may result in unreliable transmission and reception of data. The strength indicator is located at the top right of the meter's screen.

3. Message: "High" displayed

What it means: If you see High displayed on your device, your blood sugar level may be out of range (in excess of 600 mg/dL [High]).

Action to take: See "Troubleshooting" (pg. 54 ), and wash your hands. Dry your hands completely. Test again using fingertip testing. DO NOT test on your palm or forearm. If you still receive the same message, call your health care provider immediately.

4. Message: Software updates are available in the App Store for your BGM

What it means: A critical or non-critical software update is available on the App Store. If the update is critical, you won't be able to use other meter functions until the software update is completed. If the update is non-critical, you may dismiss the update. The update will appear every 24 hours and will also be available in the App Store.

Action to take: See "Apps" section (pg. 32 ). Follow the prompts to install the software update.

5. Error: Meter freezes or restarts

What it means: Your meter may require restarting or software updates.

Action to take: Check your meter's battery charge level. Charge your meter if low or very low battery. Turn off your device, wait 10 seconds, and turn your device back on. Check if your meter has pending software updates and install them. If the issue persists, contact Customer Care at (888) 666-4530 (US) or (1) (978) 610-4688 (International) (9AM-9PM ET, M-F).

Communication errors

If you are having trouble sending your results, you may want to try the following:

- If testing indoors, move closer to a window for better cellular reception
- Check to see if there is cellular coverage in your area. Two bars or less of wireless signal strength may result in unreliable transmission and reception of data.
- Restart your meter and send reading within Logs.

If your meter has had an extended period of inactivity or excessive data usage for ancillary applications, your cellular service may have been discontinued in accordance with the Terms of Use. Contact Customer Care for re-activation options.

If the above actions do not resolve your problem, contact Customer Care at (888) 666-4530 (US) or (1)(978) 610-4688 (International) (9AM-9PM ET, M-F), or e-mail VCM-customercare@philips.com.

Patient Consent

Upon initial boot up, you will be prompted to accept the Terms of Use included in this user manual. You will not be able to start a test or use the meter functions until the terms have been accepted. If you do not wish to accept the terms, contact Customer Care at (888) 666-4530 (US) or (1)(978) 610-4688 (International) (9AM-9PM ET, M-F), or e-mail VCM-customercare@philips.com for opt-out options.

Terms of Warranty

Telcare, LLC, a Philips Company ("Philips") warrants the Philips Connected Blood Glucose Meter (BGM) to be free of defects in workmanship and materials under normal use for a period of one (1) year from the date of purchase.

The liability of Philips is limited to repair or replacement. In no event shall Philips be liable for any collateral or consequential damages or loss.

Instruments subjected to misuse, abuse, neglect, exposure to water, unauthorized repair, or modification are excluded from this warranty.

This guarantee specifically excludes expendables and consumables. Lost, stolen, or irretrievable devices are also specifically excluded from this warranty.

Transfer of ownership will render this warranty null and void.

Please direct any questions or claims to the authorized dealer from whom you purchased your system, or to:

Telcare, LLC, a Philips Company
600 West Ridge Road, Suite 200
Linwood, PA 19061

Phone: (888) 666-4530 (US) or
(1) (978) 610-4688 (International)
Email: VCM-customercare@philips.com

This warranty applies only to the original purchaser of the system.

Specifications

Test Strips: Philips or BioTel Care® Test Strips

Minimum Sample: 0.8 microliters

Test Range: 20-600 mg/dL

Calibration: Plasma, calibrated against the YSI Glucose Analyser (YSI, inc.) reference instrument, which is traceable to NIST SRM 917c reference standard.

Hematocrit Range: 30-55%

Display Type: Color LCD Screen

Memory: 10,000 blood test results with date and time stamp

Result Averaging: 7, 14, 30, and 90 days averaging displayed on the BGM

Dimensions: 102L x 59.5W x 13.9H (mm)

Weight: 92 grams

Battery: Rechargeable, Lithium Ion 980mAh battery, replaceable.

Battery Life: 1 year, call Customer Care to replace your battery or if you have battery life issues.

Automatic Sleep Mode: After 1 minute of non use (2 minutes when on the Apply Blood screen)

Operating Temperature: 50°F-104°F (10°C-40°C)

Relative Humidity: Less than 85%

Storage Conditions: Meter at 39°F-122°F (4°C~50°C). Test Strips at 39°F-86°F (4°C-30°C). Control Solution at 39°F-86°F (4°C-30°C)

J. About your system

Contact Customer Care

Telcare, LLC, a Philips Company

600 West Ridge Road, Suite 200

Linwood, PA 19061

Phone: (888) 666-4530 (US) or
(1) (978) 610-4688 (International)

Hours: M - F, 9AM - 9PM ET

E-mail: VCM-customercare@philips.com

Website: www.philips.com/virtualcaremanagement

For any medical questions, contact your health care professional.

Product and safety information

Information on health

The U.S. Food and Drug Administration (FDA) and the U.S. Federal Communications Commission (FCC) published statements and questions and answers concerning mobile devices and health. Philips encourages you to visit these websites for updated information. The BGM is not a cell phone but is considered a similar wireless device. You can access the FDA website at <https://www.fda.gov/radiation-emitting-products/home-business-and-entertainment-products/cell-phones> and the FCC website at transition.fcc.gov/oet/rfsafety/rf-faqs.html. Additional health related information is available from the World Health Organization (WHO) at www.who.int/mediacentre/factsheets/fs193/en/ and The National Cancer Institute ("NCI") at <http://www.cancer.gov/publications/blogs-newsletters/cancer-bulletin/>. In the event that you are concerned about possible health effects, the FDA suggests that you limit your own or your children's radio frequency (RF) exposure by limiting the length of your calls or by using hands free devices.

Implanted Medical Devices

To avoid potential interference, manufacturers of implanted medical devices recommend a minimum separation of six (6) inches (15.3 centimeters) between a wireless device and the medical device. Persons who have such devices should:

- Always keep the wireless device more than six (6) inches (15.3 centimeters) from the medical device.
- Do not carry the wireless device in a breast pocket .
- Hold the wireless device to the ear opposite the medical device.
- Turn the wireless device off if there is any reason to suspect that interference is taking place.
- Follow the manufacturer directions for the implanted medical device.

If you have any questions about using your wireless device with an implanted medical device, consult your health care provider.

Wireless Technology Information

The BGM is a wireless medical device which encompasses wireless technologies and functions. This device is a data only device and does not incorporate voice service. The BGM meets the Verizon Wireless standards and has been certified to function on their networks. You can only use your device on the LTE network.

This EUT (Equipment Under Test) has been shown to be capable of complying with the requirements set forth in the following standards:

- FCC 47 CFR Part 2, 22(H), 24(E)
- FCC Part 15, Subpart B
- ANSI/TIA/EIA-603-C-2004
- FCC KDB 971168 D01 Power Meas. License Digital Systems v01

This device uses 4G LTE. The following are operating characteristics of the BGM:

- Effective RF Radiated Power Output and Operating Range: Less than 24 dBm

- Modulation: QPSK. (Quadrature Phase-Shift Keying)
- Antenna Type: PIFA Antenna
- Bandwidth of Receiving Section - Is in accordance with FCC Rule Part 4E for Frequency Range
- Part 24E = 1851.25-1908.75 (mHz)

Electromagnetic Emitted Interference According to IEC 61326

The BGM is suitable for operation in the indicated electromagnetic environment. The customer and/or operator of the BGM should make sure that it is used in an electromagnetic environment as described.

IEC 61326-1: 2013, IEC 61326-2-6:2013

Test Condition	Compliance
Radiated Emission (30MHz to 1000MHz)	Complies
Conducted Emission on AC 150kHz to 30MHz	Complies
Harmonic Class A	Complies
Flicker	Complies

IEC 61326-1: 2013, IEC 61326-2-6:2013

Immunity Tests	Compliance
Electrostatic Discharge (IEC 61000-4-2) $\pm 4\text{kV}$ $\pm 8\text{kV}$	Complies
Radiated Immunity (IEC 61000-4-3) 80MHz to 1000MHz 3V/m (rms) 1400MHz to 2000MHz 3V/m (rms) 2000MHz to 2700MHz 1V/m (rms)	Complies
Electrical Fast Transient (IEC 61000-4-4) $\pm 1\text{kV}$	Complies
Surge (IEC 61000-4-5) $\pm 0.5\text{kV}$ (Line to line) $\pm 1\text{kV}$ (Line to ground)	Complies
Conducted Immunity (IEC 61000-4-6) 150KHz to 80MHz 3Vrms	Complies
Voltage Dips and Interruption (IEC 61000-4-11) 0%, 70% of Ut	Complies

FCC Interference Statement

Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 18 of the FCC Rules. These limits are designed to provide reasonable protection

against harmful interference in a residential installation. This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment when used with Philips accessories supplied or designated for this product. Use of other accessories may not ensure compliance with FCC RF exposure guidelines. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- This device complies with Part 18 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
- FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

RF Exposure Information, Specific Absorption Rate (SAR)

This device meets the government's requirements for exposure to radio waves.

This device is designed and manufactured not to exceed the emission limits for exposure to radio frequency (RF) energy set by the Federal Communications Commission of the U.S. Government.

The exposure standard for wireless mobile devices employs a unit of measurement known as the Specific Absorption Rate,

or SAR. The SAR limit set by the FCC is <1.6 W/kg. Tests for SAR are conducted using standard operating positions accepted by the FCC with the device transmitting at its highest certified power level in all tested frequency bands. Although the SAR is determined at the highest certified power level, the actual SAR level of the device while operating can be well below the maximum value. This is because the device is designed to operate at multiple power levels so as to use only the power required to reach the network. In general, the closer you are to a wireless base station, the lower the power output.

The highest SAR value for the model device as reported to the FCC when tested for use when worn on the body, as described in this user guide, is 1.07 W/kg (body-worn measurements differ among device models, depending upon available accessories and FCC requirements).

While there may be differences between the SAR levels of various devices and at various positions, they all meet the government requirement.

The FCC has granted an Equipment Authorization for this model device with all reported SAR levels evaluated as in compliance with the FCC RF exposure guidelines. SAR information on this model device is on file with the FCC and can be found under the Display Grant section of www.fcc.gov/oet/ea/fccid after searching on Grantee Code: 2AVLM and Product Code: CARE01.

For body worn operation, this device has been tested and meets the FCC RF exposure guidelines for use with an accessory that contains no metal and that positions the handset a minimum of 1.5 cm from the body. Use of other accessories may not ensure compliance with FCC RF exposure guidelines. If you do not use a body-worn accessory and are not holding the device at the ear, position the handset a minimum of 1.5 cm from your body when the device is switched on.

K. Terms of Use

Last Updated August 1, 2023

THESE TERMS OF USE SHALL BE DEEMED TO BE ACCEPTED UPON YOUR OPENING OF THE PACKAGING AND USE OF THE METER.

The following terms constitute an agreement between you and Telcare Medical Supply, LLC, a Philips Company ("Philips"), its subsidiaries and affiliates (collectively, the "Company" or "Philips" or "we" or "us"), the operator of www.philips.com/virtualcareportal and www.philips.com/virtualcaremanagement, and related websites (collectively, the "Sites"). This Agreement governs your use of the Sites, both as a visitor and a registered member as described below, and purchase, use and/or access to the Philips service and products (collectively, the "Service").

By accessing or using the Sites and the Service, you agree to be bound by the terms of this Agreement. To the extent you register on the Sites as a member, then in addition to the terms that are applicable to all users of the site, the terms that apply only to members shall apply to you as well.

Changes to Agreement

Philips may update this Agreement at its sole discretion, and you will be asked to review and agree to the revised version of the Agreement once it becomes effective. If you do not agree to a revised Agreement, you may not be permitted to continue using the Sites and the Service. If at any time you are no longer able to comply with the terms of the then-current Agreement, you must terminate this Agreement and immediately stop using the Service. Philips may change, modify, suspend, or discontinue any aspect of Service at any time. Philips may also impose limits on certain features or restrict your access to parts or all of the Service without notice or liability. You have no interest, monetary or otherwise, in any feature or content contained in the Service covered in this Agreement.

Membership Eligibility

To access certain parts of the Sites and the Service, including but not limited to the e-commerce section and forum section, you must become a member. To become a member and access the area on the Sites reserved for members (the "Member Area"), the Company requires that you are either (a) a diagnosed patient of the particular community you are joining, (b) a caregiver for a patient eligible to join such community, (c) a health care professional (e.g. doctor, nurse, health researcher, etc.), (e) registered user for forum pages, (f) registered e-commerce user, or (d) a guest as authorized by a Philips employee.

The Sites and the Services Do Not Provide Medical Advice.

All of the material provided on the sites and as part of the service, such as text, treatments, dosages, outcomes, charts, patient profiles, graphics, photographs, images, advice, messages, forum postings, and any other material provided on the sites and as part of the service are informational purposes only and are not a substitute for professional medical advice or treatment. You hereby agree and acknowledge that such information shall not constitute medical advice, directly or indirectly, whatsoever. You should always seek the advice of your physician or other qualified health provider with any questions you may have regarding your health. Never disregard professional medical advice or delay in seeking it because of something you have read on the sites or received as part of the service.

If you think you may have a medical emergency, call your doctor or 911 immediately. Philips does not recommend or endorse any specific tests, physicians, products, procedures, opinions, or other information that may be mentioned in on the sites or as part of the service. Reliance on any information provided by Philips, by persons appearing on the sites or at the part of the service at the invitation of Philips, or by other members is solely at your own risk.

Acceptable and Lawful Use of the Sites by Members

Members shall not post or upload any information or other content on the Sites that (a) is false, inaccurate or misleading; (b) is obscene or indecent; (c) infringes any copyright, patent,

trademark, trade secret or other proprietary rights or rights of publicity or privacy of any party; or (d) is defamatory, libelous, threatening, abusive, hateful, or contains pornography.

Members shall not interfere with other members' use and enjoyment of Sites (for example, by harassing, threatening or stalking another participant of a chat room, forum, seminar, or threaded message, including the use of expressions of bigotry, racism, hatred, or profanity). Members may not use the Sites to conduct any activity that is illegal or violates the rights of others, provide instructional information about illegal activities, or promote physical harm or injury against any group or individual. All members represent and warrant that the information they provided when registering as a member, and all information that they subsequently provide regarding themselves and their membership, is true and accurate and not misleading.

Use of Sites by Members and Non-Members

You may not use any robot, spider, scraper, or other automated means to access the Sites, the Service or content or services provided on the Sites for any purposes. You may not post content on the Sites or as part of the Service that contains any viruses, Trojan horses, worms, time bombs, spiders, or other computer programming routines that are intended to damage, detrimentally interfere with, surreptitiously intercept or expropriate any system, data or personal information. You shall not attempt to make the Sites or the Service unavailable through denial-of-service attacks or similar means. You shall not use contact information provided by members, or collect information about our members, to facilitate the sending of unsolicited bulk communications such as SPAM or SPIM or allow others use of your membership account to take such actions.

Non-Commercial Use by Members

The member area of the Sites and the Service (the "Member Area") and the content and information contained in such Member Area is for the personal use of individual members only and may not be used in connection with any commercial endeavors. Organizations, companies, and/or businesses may not become members and should not use the Sites or any part

of the Service without the prior express written consent from the Company.

Privacy and Your Data

You agree that you have read, understood and accept the terms of Philips's Privacy Policy. This policy governs the collection, use, and sharing of personal and non-personal information from you when using the Sites and the Service. As detailed in our Notice of Privacy Practices, Philips may use your medical information for the following purposes:

- Treatment. We may use or disclose medical information about you to provide you with medical services and/or to assist medical providers in coordinating and managing your care.
- Payment. We may use or disclose your health information so that the services you receive may be billed to and payment may be collected from you, an insurance company, or a third party. We may also disclose your medical information to contractors who provide claims processing or other services to Philips.
- Health care operations. We may use or disclose your health information to perform necessary internal functions. For example, we may use your data to help us train new staff and to conduct quality improvement activities. We may combine medical information about many patients to decide what additional services to offer and whether certain services are effective.

As part of the Service, you control who may access your data, along with any authorized person who may access such data on your behalf (and grant access to such data) and other invited users. If you share your data generated by the Service, then you agree that anyone you've shared such data with may access and use that data in which such use shall not be subject to Philips's Privacy Policy. You agree that SMS/Text and e-mail messages sent to authorized users are for informational purposes only. Furthermore, by using the Service, you hereby consent to the access and review of such data by your health care insurance provider (the "Insurance Provider") and hereby waive and release the Insurance Provider and its respective affiliates from any and all damages, losses, actions, costs, and claims in connection with

the access of your data resulting from your use of the Service and in connection with the administration of your participation in the Insurance Provider's health care plan.

Furthermore, you hereby acknowledge and agree that Philips shall be permitted to view any and all of your data ("End User Data") for the sole purpose of complying with FDA regulations and its obligations therein. In addition, Philips may, in its discretion, extract from the End User Data information on an aggregated and anonymous basis, including without limitation your age, gender, and prescription protocols (the "De-Identified Information") from which identifiers have been removed pursuant to 45 C.F.R. § 164.514 (together with any other relevant statutes or regulations, as amended from time to time), such that the De-Identified Information is not Protected Health Information ("PHI") and is not subject to the privacy restrictions of HIPAA and any and all related privacy regulations. Without limiting Philips's obligation to preserve the confidentiality of your PHI, Philips may use your De-Identified Information for statistical, analytical, and any academic or commercial purposes on an aggregated basis. Philips shall not use or disclose the De-Identified Information in any manner inconsistent with HIPAA or any other statute or regulation.

You hereby agree and acknowledge that your testing data and blood glucose readings ("Testing Data") will not be automatically sent to or monitored by your physician, medical advisor or any health personnel that you seek or receive health services from in connection with your health management needs (collectively, "Health Agents"). Accordingly, you hereby waive and release any such Health Agents and any respective affiliates from any and all liability, damages, losses, actions, costs, and claims in connection with the access of your Testing Data resulting from your use of the Service.

Posting and Use of Content in the Member Area

The Member Area includes community areas, such as forums and member areas, where members may post messages, images, and other content. If you are a member, you (or the author) owns the copyright in the messages, images, and other content you post in the Member Area, but by posting such content to the Member

Area you grant Philips and our affiliates the right and irrevocable, unlimited, worldwide, and royalty-free license to use, copy, display, perform, distribute, translate, edit, and create derivative works of your postings, subject to the terms of the Privacy Policy.

While some community areas in the Member Area are monitored periodically for topicality, Philips has no obligation to prescreen postings and is not responsible for their content. We encourage you to notify us of inappropriate or illegal content and we reserve the right to remove postings for any reason.

You agree not to disclose to any person or entity personally identifiable information about other members that you learn using the Sites or the Service (whether posted in the Member Area by a member or e-mailed to you by a member) without the express consent of such member. You may disclose information of a general nature (that could not identify the member who provided such information or whom such information is about) to third parties outside the Sites and the Service, subject to the above restriction on non-commercial use.

Member Password and Login Identity

You are responsible for maintaining the confidentiality of your member password and login for the Sites and the Service and are fully responsible for all activities that occur under your password or account with or without your knowledge. If you knowingly provide your login and password information to another person, your membership privileges to the Sites and the Service may be suspended temporarily or terminated. You agree to immediately notify Philips of any unauthorized use of your membership password or login or any other breach of security.

Links to Other Websites

The Sites may contain links to other websites that are not under the control of Philips. The inclusion of any link does not imply endorsement by Philips of such website and Philips is not responsible for the content of such linked websites.

Membership Termination

You agree that the Company may, with or without cause, immediately terminate your Philips membership and access to the Member Area, the Sites and the Service without prior notice.

Without limiting the foregoing, the following may lead to a termination by the Company of a member's membership and use of the Sites and the Service: (a) breaches or violations of this Agreement or other incorporated agreements or guidelines, (b) requests by law enforcement or other government agencies, (c) a request by you (self-initiated membership cancellation), (d) unexpected technical issues or problems, and (e) extended periods of inactivity (device has not transmitted a test reading for 60 days). Philips has no obligation to maintain, store, or transfer to you information or data that you have posted on or uploaded to the Member Area.

Modifications to this Agreement

We reserve the right to modify this Agreement at any time, and without prior notice, by posting amended terms on the Sites.

We encourage you to review this Agreement periodically for any updates or changes, which will have an updated effective date.

Fees

Membership to the Sites is currently free. Philips reserves the right to from time to time change the terms of this Agreement and charge for membership and/or to provide additional optional services for a fee in connection with the Sites and/or the Service. Products may be purchased from philips.com/virtualcaremanagement.

Purchases; Other Terms and Conditions

The Company obligations, if any, with regard to the Services are governed solely by the agreements pursuant to which they are provided, and nothing on the Sites should be construed to alter such agreements.

Prices and availability of goods and services are subject to change without notice. Philips may make changes to any Services offered on the Sites, or to the applicable prices for any such products or services, at any time, without notice. The materials and content on the Sites with respect to products and services may be out of date, and Philips makes no commitment to update the materials or content on the Sites with respect to such products and services.

The Company reserves the right to revoke any stated offer and to correct any errors, inaccuracies or omissions including after an order has been submitted and whether or not the order has been confirmed and your credit card charged. If your credit card has already been charged for the purchase and your order is canceled, Philips will issue a credit to your credit card account in the amount of the charge. Individual bank policies will dictate when this amount is credited to your account. If you are not fully satisfied with your purchase, you may return it in accordance with the Company Return Policy.

Paying for Your Order

Generally, credit and debit cards are not charged until we either ship the item(s) to you or confirm availability (at which time you will be charged only for the goods we have actually shipped along with any appropriate taxes or shipping charges). However, we may pre-authorize your order amount with your credit or debit card issuer at the time you place the order, which may have an effect on your available credit line. Please contact your credit card issuer for more information.

If you ordered a special delivery item, you will be charged once a delivery time is confirmed. For digitally delivered orders, your credit or debit card will be charged at the time that you initiate the download of the product.

Validating Your Order

After you place an order using our shopping cart, we will check the information you give us for validity, by verifying your method of payment or shipping address. We reserve the right to reject any order you place with us, and/or to limit quantities on any order, without giving any reason. If we reject your order, we will attempt to notify you using the e-mail address you have given us with the order.

Your credit or debit card will normally not be charged if we reject an order, but we will process a refund if the charge has been made against your card.

Order Acceptance/Confirmation

Your receipt of an electronic or another form of order confirmation does not signify our acceptance of your order, nor

does it constitute confirmation of our offer to sell. Philips reserves the right at any time after receipt of your order to accept or decline your order for any reason.

Return Policy

All sales made through the Sites are subject to Company return policy. For a complete explanation of the return policy, please see our Return Policy.

Order Limitations

We reserve the right to reject any order you place with us, and/or to limit quantities on any order, without giving any reason. If we reject your order, we will generally attempt to notify you using the e-mail address you gave us when you placed the order.

Philips may, at its own discretion, limit or cancel quantities purchased per person, per household or per order. Philips also reserves the right to reject any order you place with us. These restrictions may include orders placed by the same Sites account, the same credit card, and orders that use the same billing and/or shipping address. In the event we make a change to an order, we will attempt to notify you by contacting the e-mail and/or billing address provided at the time the order was made. Philips reserves the right to limit or prohibit sales to distributors.

Shipping

All orders placed on the Sites are subject to product availability and will be shipped according to Company shipping policies. In the rare event that your order never arrives, notification of lost items must be received within 60 days from receipt of the shipping confirmation e-mail.

Please review the Shipping page for additional information on shipping times, rates, and policies.

Disclaimer of Warranties

The sites, their applicable the content and services are provided on an "as is" and "as available" basis. Philips does not make, and disclaims all express and implied warranties and representations, including, but not limited to, any implied warranty of fitness for a particular purpose, with regard to the sites and the service, the sites content or any advice or services provided through the sites to the extent permitted by law. Philips does not warrant

that access to the sites, their applicable content, or services will be uninterrupted or error-free or that defects in the sites or the services will be corrected.

The advice, recommendations, information, and conclusions posted or e-mailed by other members of the sites are not in any way vetted, approved or endorsed by Philips, and you use such information at your own risk.

Not all products sold on the sites have warranties. Warranties are limited to individual products when applicable. Refer to the product for warranty information.

Limitation of Liability.

Except where prohibited by law, under no circumstances shall Philips be liable for any indirect, punitive, incidental, special, or consequential damages (even if it has been advised of the possibility of such damages) due to your use of the sites, the services or due to your reliance on any of the content contained or the services.

Indemnification

You agree to indemnify and hold harmless the Company, its officers, employees, agents, subsidiaries, affiliates, and other partners, from and against any claims, actions or demands, liabilities and settlements including without limitation, reasonable legal and accounting fees, resulting from, or alleged to result from, your violation of this Agreement.

International Users

Philips makes no claims that information on the Sites is appropriate or may be used or downloaded outside of the United States. Personal information ("Information") that is submitted to this Sites and the Service will be collected, processed, stored, disclosed and disposed of in accordance with applicable U.S. law and our Privacy Policy. If you are a non-U.S. user, you acknowledge and agree that the Company may collect and use your Information and disclose it to other entities outside your resident jurisdiction. In addition, such Information may be stored on servers located outside your resident jurisdiction. U.S. law may not provide the degree of protection for Information that is available in other countries. By providing us with your

Information, you acknowledge that you consent to the transfer of such Information outside your resident jurisdiction as detailed in our Privacy Policy. If you do not consent to such transfer, you may not use these Sites or the Service.

Member Notices

If you register as a member, you agree that Philips may send notices to you by e-mail at the e-mail address you provide when registering to become a member (or which you later update using the functionality of the Sites or the Service).

Governing Law and Venue

This Agreement shall be governed by the laws of the state of Delaware and the applicable federal laws of the United States of America. All disputes arising under, or in any way connected with membership in or use of the Sites, shall be litigated exclusively in the state and federal courts residing in the state of Delaware, and in no other court or jurisdiction.

You hereby submit to the jurisdiction of the state and federal courts sitting in the state of Delaware.

Miscellaneous Terms

Philips and the Philips logo are trademarks of Koninklijke Philips N.V. You agree not to display or use these trademarks in any manner without prior, written permission from Philips.

Philips may assign this Agreement at any time to a subsidiary or parent company or to a successor to its business as part of a merger or sale of substantially all of its assets. You may not assign or transfer this Agreement.

If any provision of this Agreement is held to be unenforceable for any reason, the remaining provisions will be unaffected and remain in full force and effect.

Some of the treatment information provided on the Sites and the Service is from third party sources. Every effort has been made to ensure that the information provided is accurate, up-to-date, and complete, but no guarantee is made to that effect. The absence of a warning for a given drug or drug combination in no way should be construed to indicate that the drug or drug combination is safe, effective or appropriate for any given patient.

Telcare, LLC, a Philips Company
600 West Ridge Road, Suite 200
Linwood, PA 19061, USA

**Phone: (888) 666-4530 (US) or
(1) (978) 610-4688 (International)**

VCM-customercare@philips.com

www.philips.com/virtualcaremanagement

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© 2023 Koninklijke Philips N.V. All rights reserved.
**Telcare, LLC, una empresa de Philips