


# LOOPING EVENT

Event Monitors



## SPECIFICATIONS

- 5 Event Memory
- 30 Day Battery Life
- Record 60 pre/30 post
- Lead II ECG
- Lead Loss Detection
- Mute Function
- Small/Compact Design
- Memory Full Indicator
- 4X Transmission Speed

The BioTel Heart Looping Event is a patient-activated loop event monitor that offers the technology to document symptomatic cardiac events.

The simple single button operation allows patients to document and record transient symptoms including fatigue, palpitations, dizziness, chest discomfort and shortness of breath.