

# CHEST PLATE

Post Symptom Event Monitor



The BioTel Heart Chest Plate is a non-looping cardiac event monitor that provides the technology to document post-symptomatic events.

Patients experiencing sustained symptomatic events are afforded ease of use with the small compact design of the chest plate. The chest plate records a 60-second event when pressed against the chest.

## CHEST PLATE SPECIFICATIONS

- Transmits via Landline
- No Lead Wires
- Easy online enrollment with every patient activated event reported
- Compact Design
- Symptomatic Patients Only
- 30-day Battery Life



CardioNet & LifeWatch