

DUAL ALERT AFIB

Auto Detection Loop Event Monitor



DUAL ALERT AFIB SPECIFICATIONS

- Field Programmable
- 6 Event Memory
- 30 Day Battery Life
- Record 60 pre/30 post
- Lead II ECG
- Two or Four Lead Hook-Up
- Bradycardia 30 to 60 bpm
- Tachycardia 120 to 220 bpm
- Cardiac Pause 2.0 to 3.5 secs
- Atrial Fibrillation Onset

The BioTel Heart Dual Alert AFIB is an auto-detection cardiac event monitor providing capture of both symptomatic and asymptomatic cardiac events.

An integrated auto detection algorithm allows the Dual Alert AFIB to trigger with or without patient activation on bradycardia, tachycardia, pause and onset of atrial fibrillation.