

PHILIPS

BioTel
HEART



Delivering necessary cardiac monitoring – without an office visit

Home Connect™ Program: A telehealth initiative to support direct to patient arrhythmia monitoring

Philips BioTel Heart's Home Connect solution is ideally suited to support your practice with a proven at-home patient activation model.

We have over 20 years of experience with our Home Connect model and a robust logistics network to ensure that patients receive their at-home activation kit quickly.

Direct to patient arrhythmia monitoring



Deliver access to care remotely



Simple set up for your staff and patients



View and interpret reports anywhere you have an internet connection

Simple patient set-up



Step-by-step patient education guide included in every device kit



24/7 access to online instructional video to assist patients



Patient Services team available to support patients

“The workflow is seamless and efficient for our staff and enables continuity of care for our patients in a safe manner. Due to the Home Connect program, our practice easily transitioned to telehealth which will remain a core component of how we deliver care now and in the future.”


Michelle, Exercise Physiologist – A large cardiology practice serving both urban and rural areas

Home Connect solutions

MCOT® (Mobile Cardiac Outpatient Telemetry)

Up to 30 days

- SmartDetectAI algorithm delivers near-real-time arrhythmia detection with rate, rhythm, AF with p-wave analysis and QRS morphology
- Proven to detect AF with 100% sensitivity and 100% positive predictivity in the detection of ≥ 30 -second AF episodes¹—no false AF negatives or positives¹
- Superior detection and high diagnostic yield expedite patient treatment decisions
- 2 channels of ECG data for a more complete picture
- Multiple wear options enable patients to switch to an alternate form factor for continuous monitoring.




Patient wears near real-time MCOT for up to 30 days. Notifications sent per provider protocol.[†]

ePatch® extended wear Holter

Up to 14 days

- Easy-to-apply devices your patients are willing to wear for as long as you prescribe — up to 14 days²
- Comprehensive reports driven by advanced AI software that captures vital cardiac data to make informed care decisions
- Reports delivered in less than 48 hours



Patient wears extended Holter for up to 14 days.

Learn more about our Home Connect offerings or begin enrollment by contacting your local Account Executive, calling 866-426-4401, or visiting gobio.com/telehealth

1000 Cedar Hollow Road, Malvern, PA 19355
www.philips.com • 866-426-4402

LifeWatch, CardioNet and BioTel Heart are trademarks of BioTelemetry, Inc. BioTelemetry reserves the right to change specifications at any time without notice. © 2022 BioTelemetry. All rights reserved. Doc. 220-0864-01 Rev. B



1. Based on MIT-BIH (Massachusetts Institute of Technology- Beth Israel Hospital) Arrhythmia Database testing of ≥ 30 -second AF episodes. (FDA 510k submission)
2. Patch must be changed within 5 days or sooner.
† Per provider predetermined protocol.
© 2022 Koninklijke Philips N.V. All rights reserved.